

Human SPTAN1 Knockdown Cell Line (WB-Validated)



Catalog #: C64875

Aliases

SPTAN1; Spectrin Alpha, Non-Erythrocytic 1; Spectrin Alpha Chain, Non-Erythrocytic 1; Spectrin, Non-Erythroid Alpha Subunit; Fodrin Alpha Chain; Alpha-II Spectrin; Alpha-Fodrin; SPTA2; Epididymis Secretory Sperm Binding Protein; Spectrin, Non-Erythroid Alpha Chain; HMND11; DEVEP; EIEE5; HMN11; SPG91; DEE5; NEAS

Background

Gene Name: SPTAN1
NCBI Gene Entry: [6709](#)

Storage

Store at liquid nitrogen for 1 year.

Kit Components

1. Human SPTAN1 Knockdown Cell Line (Wb-Validated)
2. Wild-type cell line

Parental Cell Line

Human cell line supplied by the client

Validation Methods

RT-qPCR, Western blotting (WB)

Shipping

Shipped on Dry Ice.

Instructions For Use

This knockdown cell line should be paired with wild-type cell line for use.

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

RT-qPCR analysis. HeLa cells were infected with SPTAN1-specific shRNA lentiviral particles, total RNA was extracted from wild-type and knockdown cells, RT-qPCR was performed using gene-specific primers. $\Delta Ct (Ct_{KD} - Ct_{WT})$ was used to calculate mRNA reduction (%) between wild-type and knockdown cells using the following formula: $(1 - 1/2^{\Delta Ct}) \times 100\%$.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis. SPTAN1 protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting. Hsp90 α served as a loading control. The blots were incubated with primary antibodies against SPTAN1 and Hsp90 α, respectively, followed by incubating with HRP-conjugated goat anti-mouse secondary antibody. Images were developed using FeQ™ ECL Substrate Kit.