

Anti-CD1a Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0001-APC

Aliases

CD1A; CD1a Molecule; T-Cell Surface Glycoprotein CD1a; T-Cell Surface Antigen T6/Leu-6; CD1A Antigen, A Polypeptide; HTa1 Thymocyte Antigen; CD1a Antigen; CD1; Epidermal Dendritic Cell Marker CD1a; Differentiation Antigen CD1-Alpha-3; Cortical Thymocyte Antigen CD1A; Cluster Of Differentiation 1 A; FCB6; HTA1; R4; T6

Background

Gene Name: CD1A

NCBI Gene Entry: [909](#)

UniProt Entry: [P06126](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14395

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD1a Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM