

Anti-CD4 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0005-APC

Aliases

CD4; CD4 Molecule; T-Cell Surface Glycoprotein CD4; Leu-3; T4; T-Cell Surface Antigen T4/Leu-3; CD4 Antigen (P55); CD4 Receptor; CD4 Antigen; CD4mut; IMD79; OKT4D

Background

Gene Name: CD4

NCBI Gene Entry: [920](#)

UniProt Entry: [P01730](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14415

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD4 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM