

# Anti-CD31 Mouse Monoclonal Antibody, APC Conjugated



**Catalog #: F0040-APC**

## Aliases

PECAM1; Platelet And Endothelial Cell Adhesion Molecule 1; CD31 Antigen; CD31; Platelet Endothelial Cell Adhesion Molecule; PECAM-1; EndoCAM; GPIIA'; PECA1; Platelet/Endothelial Cell Adhesion Molecule 1; Platelet Endothelial Cell Adhesion Molecule-1; CD31/EndoCAM

## Background

Gene Name: PECAM1

NCBI Gene Entry: [5175](#)

UniProt Entry: [P16284](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14590

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD31 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)