

Anti-CD43 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0050-APC

Aliases

SPN; Sialophorin; GPL115; CD43; Leukosialin; LEU-22; LSN; Sialophorin (GpL115, Leukosialin, CD43); PLeukocyte Sialoglycoprotein; Galactoglycoprotein; GALGP; CD43 Antigen

Background

Gene Name: SPN

NCBI Gene Entry: [6693](#)

UniProt Entry: [P16150](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14640

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD43 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM