

Anti-CD50 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0058-APC

Aliases

ICAM3; Intercellular Adhesion Molecule 3; ICAM-R; CDW50; CD50; ICAM-3; CD50 Antigen; CDw50

Background

Gene Name: ICAM3

NCBI Gene Entry: [3385](#)

UniProt Entry: [P32942](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14680

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD50 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM