

Anti-CD52 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0059-APC

Aliases

CD52; CD52 Molecule; HE5; EDDM5; CDW52; Human Epididymis-Specific Protein 5; CD52 Antigen (CAMPATH-1 Antigen); Epididymal Secretory Protein E5; Cambridge Pathology 1 Antigen; CAMPATH-1 Antigen; Epididymis Secretory Sperm Binding Protein Li 171mP; CDW52 Antigen (CAMPATH-1 Antigen); CD52 Antigen; HEL-S-171mP; CDw52; He5

Background

Gene Name: CD52

NCBI Gene Entry: [1043](#)

UniProt Entry: [P31358](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14685

Species Reactivity: Human, cynomolgus, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD52 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM