

Anti-CD66e Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0070-APC

Aliases

CEACAM5; CEA Cell Adhesion Molecule 5; Carcinoembryonic Antigen-Related Cell Adhesion Molecule 5; CD66e; CEA; Carcinoembryonic Antigen Related Cell Adhesion Molecule 5; Meconium Antigen 100; Carcinoembryonic Antigen; CD66e Antigen

Background

Gene Name: CEACAM5

NCBI Gene Entry: [1048](#)

UniProt Entry: [P06731](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14740

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD66e Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM