

# Anti-CCD99 Mouse Monoclonal Antibody, APC Conjugated



**Catalog #: F0082-APC**

## Aliases

CD99; CD99 Molecule (Xg Blood Group); CD99 Antigen; MIC2X; MIC2Y; MIC2; Antigen Identified By Monoclonal Antibodies 12E7, F21 And O13; T-Cell Surface Glycoprotein E2; E2 Antigen; Antigen Identified By Monoclonal 12E7, Y Homolog; MIC2 (Monoclonal Antibody 12E7); Cell Surface Antigen HBA-71; Cell Surface Antigen 12E7; Cell Surface Antigen O13; Surface Antigen MIC2; CD99 Molecule; Protein MIC2; HBA71; MSK5X; 12E7

## Background

Gene Name: CD99

NCBI Gene Entry: [4267](#)

UniProt Entry: [P14209](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14800

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG2a

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CCD99 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)