

Anti-CD105 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0085-APC

Aliases

ENG; Endoglin; END; HHT1; CD105 Antigen; CD105; ORW1; Osler-Rendu-Weber Syndrome 1; Soluble Endoglin; ORW

Background

Gene Name: LAMP1

NCBI Gene Entry: [2022](#)

UniProt Entry: [P17813](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14815

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD105 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM