

Anti-CD106 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0086-APC

Aliases

VCAM1; Vascular Cell Adhesion Molecule 1; CD106; Vascular Cell Adhesion Protein 1; CD106 Antigen; INCAM-100; V-CAM 1; VCAM-1

Background

Gene Name: VCAM1

NCBI Gene Entry: [7412](#)

UniProt Entry: [P19320](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14820

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD106 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM