

# Anti-CD115 Rat Monoclonal Antibody, APC Conjugated



**Catalog #: F0087-APC**

## Aliases

CSF1R; Colony Stimulating Factor 1 Receptor; C-FMS; CD115; CSFR; FMS; McDonough Feline Sarcoma Viral (V-Fms) Oncogene Homolog; Macrophage Colony-Stimulating Factor 1 Receptor; Proto-Oncogene C-Fms; CSF-1 Receptor; CD115 Antigen; EC 2.7.10.1; M-CSF-R; CSF-1R; Macrophage Colony Stimulating Factor I Receptor; FMS Proto-Oncogene; EC 2.7.10; BANDDOS; CSF-1-R; HDLS1; FIM2; GPSC; HDLS

## Background

Gene Name: CSF1R

NCBI Gene Entry: [1436](#)

UniProt Entry: [P07333](#)

## Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB14825

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

## Immunogen

The immunogen information is proprietary.

## Isotype

Rat IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD115 Rat Monoclonal Antibody, APC Conjugated

PAGE 2

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)