

Anti-CD140b Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0095-APC

Aliases

PDGFRB; Platelet Derived Growth Factor Receptor Beta; PDGFR1; CD140b; JTK12; PDGFR; Platelet-Derived Growth Factor Receptor, Beta Polypeptide; Beta-Type Platelet-Derived Growth Factor Receptor; Platelet-Derived Growth Factor Receptor Beta; Platelet-Derived Growth Factor Receptor 1; CD140 Antigen-Like Family Member B; PDGF-R-Beta; EC 2.7.10.1; PDGFR-Beta; PDGFR-1; Beta Platelet-Derived Growth Factor Receptor; Activated Tyrosine Kinase PDGFRB; CD140b Antigen; NDEL1-PDGFRB; EC 2.7.10; IBGC4; PENTT; IMF1; KOGS

Background

Gene Name: PDGFRB

NCBI Gene Entry: [5159](#)

UniProt Entry: [P09619](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14865

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD140b Mouse Monoclonal Antibody, APC Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM