

Anti-CD180 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0101-APC

Aliases

CD180; CD180 Molecule; RP105; CD180 Antigen; Ly78; LY64; Radioprotective 105 KDa Protein; Lymphocyte Antigen 64; Lymphocyte Antigen 64 (Mouse) Homolog, Radioprotective, 105kD; Lymphocyte Antigen-64, Radioprotective, 105kDa

Background

Gene Name: CD180

NCBI Gene Entry: [4064](#)

UniProt Entry: [Q99467](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14895

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD180 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM