

Anti-CD274 Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0108-APC

Aliases

CD274; CD274 Molecule; PD-L1; B7H1; PDL1; Programmed Cell Death 1 Ligand 1; CD274 Antigen; B7 Homolog 1; PDCD1LG1; B7-H1; B7-H; PDCD1 Ligand 1; PDCD1L1; HPD-L1; Programmed Death Ligand 1

Background

Gene Name: CD274

NCBI Gene Entry: [29126](#)

UniProt Entry: [Q9NZQ7](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14930

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD274 Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM