

Anti-CD307e Mouse Monoclonal Antibody, APC Conjugated



Catalog #: F0112-APC

Aliases

FCRL5; Fc Receptor Like 5; BXMAS1; FCRH5; IRTA2; CD307e; Immune Receptor Translocation-Associated Protein 2; Fc Receptor-Like Protein 5; Fc Receptor Homolog 5; FcR-Like Protein 5; Immunoglobulin Superfamily Receptor Translocation Associated 2 (IRTA2); Fc Receptor-Like 5; CD307e Antigen; PRO820; CD307; FcRL5; FcRH5

Background

Gene Name: FCRL5

NCBI Gene Entry: [83416](#)

UniProt Entry: [Q96RD9](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14950

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD307e Mouse Monoclonal Antibody, APC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM