

Anti-CD22 Mouse Monoclonal Antibody, Biotin Conjugated



Catalog #: F0134-B

Aliases

CD22; CD22 Molecule; SIGLEC2; CD22 Antigen; SIGLEC-2; Sialic Acid-Binding Ig-Like Lectin 2; B-Lymphocyte Cell Adhesion Molecule; T-Cell Surface Antigen Leu-14; B-Cell Receptor CD22; BL-CAM; Sialic Acid Binding Ig-Like Lectin 2; Siglec-2

Background

Gene Name: CD22

NCBI Gene Entry: [933](#)

UniProt Entry: [P20273](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14535

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Biotin

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

0.25-4 µg/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD22 Mouse Monoclonal Antibody, Biotin Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM