

Anti-CD326 Mouse Monoclonal Antibody, Biotin Conjugated



Catalog #: F0145-B

Aliases

EPCAM; Epithelial Cell Adhesion Molecule; TROP1; KSA; Tumor-Associated Calcium Signal Transducer 1; GA733-2; Ber-Ep4; TACSTD1; Ep-CAM; MOC-31; BerEp4; EGP40; EGP-2; KS1/4; MIC18; MK-1; M4S1; ESA; Major Gastrointestinal Tumor-Associated Protein GA733-2; Trophoblast Cell Surface Antigen 1; Adenocarcinoma-Associated Antigen; Cell Surface Glycoprotein Trop-1; Epithelial Glycoprotein 314; TACST-1; HEA125; 323/A3; CO-17A; EGP314 ; EGP34; CD326; MOC31; 17-1A; Ly74; MH99; Membrane Component, Chromosome 4, Surface Marker (35kD Glycoprotein); Antigen Identified By Monoclonal Antibody AUA1; Human Epithelial Glycoprotein-2; Epithelial Cell Surface Antigen; Epithelial Glycoprotein; KS 1/4 Antigen; CD326 Antigen; HEGP314; HNPCC8; LYNCH8; DIAR5; M1S2; EGP

Background

Gene Name: EPCAM

NCBI Gene Entry: [4072](#)

UniProt Entry: [P16422](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14960

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: Biotin

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM


Anti-CD326 Mouse Monoclonal Antibody, Biotin Conjugated

PAGE 2

0.25-4 µg/Tests

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM