

# Anti-CD8 Mouse Monoclonal Antibody, PE-Cy5 Conjugated



**Catalog #: F0152-PC5**

## Aliases

CD8A; CD8 Subunit Alpha; T-Cell Surface Glycoprotein CD8 Alpha Chain; CD8alpha; P32; T-Lymphocyte Differentiation Antigen T8/Leu-2; CD8 Antigen, Alpha Polypeptide (P32); CD8a Molecule; CD8; Leu2 T-Lymphocyte Antigen; OKT8 T-Cell Antigen; T-Cell Antigen Leu2; T Cell Co-Receptor; T8 T-Cell Antigen; CD8a Antigen; IMD116; Leu2; MAL

## Background

Gene Name: CD8A

NCBI Gene Entry: [925](#)

UniProt Entry: [P01732](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14440

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD8 Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)