

Anti-CD11a Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0158-PC5

Aliases

ITGAL; Integrin Subunit Alpha L; LFA-1; LFA1A; CD11A; Integrin, Alpha L (Antigen CD11A (P180), Lymphocyte Function-Associated Antigen 1; Alpha Polypeptide); Lymphocyte Function-Associated Antigen 1, Alpha Polypeptide; Leukocyte Function-Associated Molecule 1 Alpha Chain; Leukocyte Adhesion Glycoprotein LFA-1 Alpha Chain; CD11 Antigen-Like Family Member A; Integrin Alpha-L; LFA-1A; Antigen CD11A (P180), Lymphocyte Function-Associated Antigen 1, Alpha Polypeptide; Integrin Gene Promoter; Antigen CD11A (P180); CD11a Antigen; LFA-1 Alpha

Background

Gene Name: ITGAL

NCBI Gene Entry: [3683](#)

UniProt Entry: [P20701](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15035

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD11a Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM