

Anti-CD25 Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0167-PC5

Aliases

IL2RA; Interleukin 2 Receptor Subunit Alpha; CD25; Interleukin-2 Receptor Subunit Alpha; Interleukin 2 Receptor, Alpha; IL-2 Receptor Subunit Alpha; IL-2R Subunit Alpha; TAC Antigen; IDDM10; IL2R; P55; Insulin-Dependent Diabetes Mellitus 10; CD25 Antigen; IL-2-RA; IL2-RA; IMD41; TCGFR

Background

Gene Name: IL2RA

NCBI Gene Entry: [3559](#)

UniProt Entry: [P01589](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14560

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD25 Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM