

Anti-CD31 Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0168-PC5

Aliases

PECAM1; Platelet And Endothelial Cell Adhesion Molecule 1; CD31 Antigen; CD31; Platelet Endothelial Cell Adhesion Molecule; PECAM-1; EndoCAM; GPIIA'; PECA1; Platelet/Endothelial Cell Adhesion Molecule 1; Platelet Endothelial Cell Adhesion Molecule-1; CD31/EndoCAM

Background

Gene Name: PECAM1

NCBI Gene Entry: [5175](#)

UniProt Entry: [P16284](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15050

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD31 Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM