

Anti-CD55 Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0177-PC5

Aliases

CD55; CD55 Molecule (Cromer Blood Group); CR; CROM; DAF; TC; CD55 Molecule, Decay Accelerating Factor For Complement (Cromer Blood Group); Complement Decay-Accelerating Factor; CD55 Antigen; Decay Accelerating Factor For Complement (CD55, Cromer Blood Group System); Cromer Blood Group Antigen; Rh Blood Group D Antigen; CHAPLE

Background

Gene Name: CD55

NCBI Gene Entry: [1604](#)

UniProt Entry: [P08174](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15075

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD55 Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM