

Anti-CD62E/CD62P Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0180-PC5

Aliases

SELE; Selectin E; CD62E; ELAM1; ESEL; Endothelial Leukocyte Adhesion Molecule 1; CD62 Antigen-Like Family Member E; Endothelial Adhesion Molecule 1; E-Selectin; LECAM2; ELAM-1; ELAM; Leukocyte Endothelial Cell Adhesion Molecule 2; Leukocyte-Endothelial Cell Adhesion Molecule 2; CD62E Antigen; Selectin-E; SELP; Selectin P; PADGEM; GMP140; CD62P; CD62; PSEL; GRMP; Platelet Activation Dependent Granule-External Membrane Protein; Selectin P (Granule Membrane Protein 140kDa, Antigen CD62); Leukocyte-Endothelial Cell Adhesion Molecule 3; CD62 Antigen-Like Family Member P; Granule Membrane Protein 140kDa; Granule Membrane Protein 140; P-Selectin; GMP-140; LECAM3; Selectin P (Granule Membrane Protein 140kD, Antigen CD62); Platelet Alpha-Granule Membrane Protein; Granulocyte Membrane Protein; CD62P Antigen; Antigen CD62; GMRP

Background

Gene Name: SELE/SELP

NCBI Gene Entry: [6401/6403](#)

UniProt Entry: [P16581/P16109](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14725

Species Reactivity: Human, pig

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD62E/CD62P Mouse Monoclonal Antibody, PE-Cy5 Conjugated


PAGE 2

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM