

Anti-CD63 Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0181-PC5

Aliases

CD63; CD63 Molecule; TSPAN30; Pltgp40; HOP-26; ME491; MLA1; AD1; Ocular Melanoma-Associated Antigen; CD63 Antigen (Melanoma 1 Antigen); Melanoma-Associated Antigen ME491; Tetraspanin-30; Granulophysin; CD63 Antigen; AD1 Antigen; Tspan-30; OMA81H; Limp1; Lysosomal-Associated Membrane Protein 3; Lysosome Integral Membrane Protein 1; Melanoma-Associated Antigen MLA1; LAMP-3

Background

Gene Name: CD63

NCBI Gene Entry: [967](#)

UniProt Entry: [P08962](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15085

Species Reactivity: Human, baboon, cynomolgus, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD63 Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM