

Anti-CD107a Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0186-PC5

Aliases

LAMP1; Lysosomal Associated Membrane Protein 1; CD107a; Lysosome-Associated Membrane Glycoprotein 1; Lysosome-Associated Membrane Protein 1; CD107 Antigen-Like Family Member A; Lysosomal-Associated Membrane Protein 1; CD107a Antigen; LGP120; LAMP-1; LAMPA

Background

Gene Name: LAMP1

NCBI Gene Entry: [3916](#)

UniProt Entry: [P11279](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15100

Species Reactivity: Human, african green, baboon, chimpanzee, cynomolgus, capuchin monkey, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD107a Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM