

Anti-TCR α/β Mouse Monoclonal Antibody, PE-Cy5 Conjugated



Catalog #: F0191-PC5

Aliases

TRA; T Cell Receptor Alpha Locus; M1-Specific T Cell Receptor Alpha Chain; TR Alpha Chain TRAV27*01J42*01C*01; TCRA, TRA@; T-Cell Antigen Receptor, Alpha Polypeptide; TR Alpha Chain TRAV38-2DV8*01J31*01C*01; T Cell Receptor Alpha Chain MC.7.G5; T-Cell Receptor, Alpha (V,D,J,C); MC.7.G5 TRA; IMD7; TCRA; TRA@; TRB; T Cell Receptor Beta Locus; T-Cell Antigen Receptor, Beta Polypeptide, T-Cell Receptor, Beta Cluster; M1-Specific T Cell Receptor Beta Chain; TR Beta Chain TRBV19*01J2S7*01C*02; TCRB, TRB@; TR Beta Chain TRBV25-1*01J2S3*01C2*01; T Cell Receptor Beta Chain MC.7.G5; T-Cell Receptor, Beta Cluster; MC.7.G5 TRB; TCRB; TRB@

Background

Gene Name: TRA/TRB

NCBI Gene Entry: [6955/6957](#)

UniProt Entry: [P0DSE1/P0DSE2](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14995

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-TCR α/β Mouse Monoclonal Antibody, PE-Cy5 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM