

Anti-CD19 Mouse Monoclonal Antibody, APC-Cy5.5 Conjugated



Catalog #: F0192-AC5.5

Aliases

CD19; CD19 Molecule; B-Lymphocyte Surface Antigen B4; T-Cell Surface Antigen Leu-12; Differentiation Antigen CD19; B-Lymphocyte Antigen CD19; CD19 Antigen; CVID3; B4

Background

Gene Name: CD19

NCBI Gene Entry: [930](#)

UniProt Entry: [P15391](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15045

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy5.5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD19 Mouse Monoclonal Antibody, APC-Cy5.5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM