

# Anti-CD20 Mouse Monoclonal Antibody, FITC Conjugated



## Catalog #: F0224-FITC

### Aliases

MS4A1; Membrane Spanning 4-Domains A1; Bp35; FMC7; CD20; B1; Membrane-Spanning 4-Domains, Subfamily A, Member 1; Leukocyte Surface Antigen Leu-16; B-Lymphocyte Antigen CD20; CD20 Antigen; Membrane-Spanning 4-Domains Subfamily A Member 1; B-Lymphocyte Cell-Surface Antigen B1; B-Lymphocyte Surface Antigen B1; CD20 Receptor; LEU-16; CVID5; S7

### Background

Gene Name: MS4A1  
NCBI Gene Entry: [931](#)  
UniProt Entry: [P11836](#)

### Application Information

Clonality: Mouse monoclonal antibody  
Clone ID: 24GB14520  
Species Reactivity: Human  
Applications Tested: Flow Cytometry (FCM)  
Conjugation: FITC

### Immunogen

The immunogen information is proprietary.

### Isotype

Mouse IgG2a

### Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

### Storage

Store at 4 °C and protected from prolonged exposure to light.

### Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM

# Anti-CD20 Mouse Monoclonal Antibody, FITC Conjugated

## Validation Data

