

Anti-CD27 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0232-FITC

Aliases

CD27; CD27 Molecule; CD27 Antigen; TNFRSF7; S152; Tp55; Tumor Necrosis Factor Receptor Superfamily, Member 7; T-Cell Activation Antigen CD27; CD27L Receptor; T14; Tumor Necrosis Factor Receptor Superfamily Member 7; T-Cell Activation Antigen S152; T Cell Activation Antigen S152; S152. LPFS2

Background

Gene Name: CD27

NCBI Gene Entry: [939](#)

UniProt Entry: [P26842](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15135

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD27 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM