

# Anti-CD27 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0233-FITC**

## Aliases

CD27; CD27 Molecule; CD27 Antigen; TNFRSF7; S152; Tp55; Tumor Necrosis Factor Receptor Superfamily, Member 7; T-Cell Activation Antigen CD27; CD27L Receptor; T14; Tumor Necrosis Factor Receptor Superfamily Member 7; T-Cell Activation Antigen S152; T Cell Activation Antigen S152; S152. LPFS2

## Background

Gene Name: CD27

NCBI Gene Entry: [939](#)

UniProt Entry: [P26842](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15140

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG2a

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

