

# Anti-CD33 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0238-FITC**

## Aliases

CD33; CD33 Molecule; SIGLEC3; CD33rSiglec; SIGLEC-3; P67; Sialic Acid-Binding Ig-Like Lectin 3; Myeloid Cell Surface Antigen CD33; CD33 Antigen (Gp67); FLJ00391; Gp67; Sialic Acid Binding Ig-Like Lectin 3; CD33 Molecule Transcript; CD33 Antigen; Siglec-3

## Background

Gene Name: CD33

NCBI Gene Entry: [945](#)

UniProt Entry: [P20138](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14600

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD33 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)