

Anti-CD35 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0241-FITC

Aliases

CR1; Complement C3b/C4b Receptor 1 (Knops Blood Group); CD35; KN; Complement Component (3b/4b) Receptor 1 (Knops Blood Group); Complement Receptor Type 1; C3b/C4b Receptor; CD35 Antigen; C3BR; Complement Component (3b/4b) Receptor 1, Including Knops Blood Group System; Complement Component 3b/4b Receptor 1 (Knops Blood Group); Knops Blood Group Antigen; Complement Receptor 1; C3-Binding Protein; C4BR

Background

Gene Name: CR1

NCBI Gene Entry: [1378](#)

UniProt Entry: [P17927](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15150

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD35 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM