

Anti-CD37 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0245-FITC

Aliases

CD37; CD37 Molecule; TSPAN26; Leukocyte Antigen CD37; Tetraspanin-26; CD37 Antigen; Tspan-26; Cell Differentiation Antigen 37; Leukocyte Surface Antigen CD37; GP52-40

Background

Gene Name: CD37

NCBI Gene Entry: [951](#)

UniProt Entry: [P11049](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15170

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD37 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM