

# Anti-CD38 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0247-FITC**

## Aliases

CD38; CD38 Molecule; ADP-Ribosyl Cyclase 1; CADPR1; ADP-Ribosyl Cyclase/Cyclic ADP-Ribose Hydrolase 1; 2'-Phospho-Cyclic-ADP-Ribose Transferase; 2'-Phospho-ADP-Ribosyl Cyclase; Cyclic ADP-Ribose Hydrolase 1; NAD(+) Nucleosidase; CD38 Antigen (P45); ADPRC1; 2'-Phospho-ADP-Ribosyl Cyclase/2'-Phospho-Cyclic-ADP-Ribose Transferase; Ecto-Nicotinamide Adenine Dinucleotide Glycohydrolase; Cluster Of Differentiation 38; CADPR Hydrolase 1; CD38 Antigen; EC 2.4.99.20; EC 3.2.2.-; EC 3.2.2.6; ADPRC1; T10

## Background

Gene Name: CD38

NCBI Gene Entry: [952](#)

UniProt Entry: [P28907](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15175

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD38 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)