

# Anti-CD41a Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0250-FITC**

## Aliases

ITGA2B; Integrin Subunit Alpha 2b; GPIIb; Integrin Alpha-IIb; PPP1R93; CD41B; CD41; GP2B; Integrin, Alpha 2b (Platelet Glycoprotein IIb Of IIb/IIIa Complex, Antigen CD41); Platelet Glycoprotein IIb Of IIb/IIIa Complex; Protein Phosphatase 1, Regulatory Subunit 93; Platelet Membrane Glycoprotein IIb; GPalpha IIb; Integrin, Alpha 2b (Platelet Glycoprotein IIb Of IIb/IIIa Complex, Antigen CD41B); Platelet Fibrinogen Receptor, Alpha Subunit; Platelet-Specific Antigen BAK; AlphaIIb Protein; CD41 Antigen; BDPLT16; BDPLT2; ITGAB; HPA3; GT1; GTA; GT

## Background

Gene Name: ITGA2B

NCBI Gene Entry: [3674](#)

UniProt Entry: [P08514](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14625

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD41a Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)