

Anti-CD42b Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0252-FITC

Aliases

GP1BA; Glycoprotein Ib Platelet Subunit Alpha; Platelet Glycoprotein Ib Alpha Chain; GPIbalpha; CD42b; Glycoprotein Ib (Platelet), Alpha Polypeptide; Antigen CD42b-Alpha; GP-Ib Alpha; GPIbA; GP1B; Platelet Membrane Glycoprotein 1b-Alpha Subunit; Mutant Platelet Membrane Glycoprotein Ib-Alpha; Platelet Membrane Glycoprotein Ib-Alpha; Glycoprotein Ib Platelet Alpha Subunit; Glycoprotein Ibalpha; CD42b Antigen; CD42b-Alpha; GPIb-Alpha; BDPLT1; BDPLT3; DBPLT3; VWDP; BSS

Background

Gene Name: GP1BA

NCBI Gene Entry: [2811](#)

UniProt Entry: [P07359](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14635

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD42b Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM