

# Anti-CD43 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0254-FITC**

## Aliases

SPN; Sialophorin; GPL115; CD43; Leukosialin; LEU-22; LSN; Sialophorin (GpL115, Leukosialin, CD43); PLeukocyte Sialoglycoprotein; Galactoglycoprotein; GALGP; CD43 Antigen

## Background

Gene Name: SPN

NCBI Gene Entry: [6693](#)

UniProt Entry: [P16150](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15180

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD43 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)