

Anti-CD45RA Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0260-FITC

Aliases

PTPRC; Protein Tyrosine Phosphatase Receptor Type C; T200; GP180; CD45R; B220; CD45; LCA; LY5; Receptor-Type Tyrosine-Protein Phosphatase C; Glycoprotein 180; CD45 Antigen; LCA; Protein Tyrosine Phosphatase, Receptor Type, C Polypeptide; T200 Leukocyte Common Antigen; Leukocyte-Common Antigen; Leukocyte Common Antigen; T200 Glycoprotein; EC 3.1.1.3.48; IMD105

Background

Gene Name: PTPRC
NCBI Gene Entry: [5788](#)
UniProt Entry: [P08575](#)

Application Information

Clonality: Mouse monoclonal antibody
Clone ID: 24GB14655
Species Reactivity: Human, chimpanzee
Applications Tested: Flow Cytometry (FCM)
Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD45RA Mouse Monoclonal Antibody, FITC Conjugated

Validation Data

