

Anti-CD47 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0263-FITC

Aliases

CD47; CD47 Molecule; IAP; Antigenic Surface Determinant Protein OA3; Leukocyte Surface Antigen CD47; MER6; OA3; CD47 Antigen (Rh-Related Antigen, Integrin-Associated Signal Transducer); Antigen Identified By Monoclonal Antibody 1D8; Integrin Associated Protein; Rh-Related Antigen; CD47 Glycoprotein; Integrin-Associated Signal Transducer; Integrin-Associated Protein; Protein MER6; CD47 Antigen

Background

Gene Name: CD47

NCBI Gene Entry: [961](#)

UniProt Entry: [Q08722](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14665

Species Reactivity: Human, baboon, cynomolgus, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD47 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM