

Anti-CD49f Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0269-FITC

Aliases

ITGA6 □ Integrin Subunit Alpha 6 □ VLA-6 □ ITGA6A □ ITGA6B □ CD49f □ CD49 Antigen-Like Family Member F □ Very Late Activation Protein 6 □ Integrin, Alpha 6 □ Integrin Alpha-6 □ Integrin Alpha6B □ CD49f Antigen □ JEB6

Background

Gene Name: ITGA6

NCBI Gene Entry: [3655](#)

UniProt Entry: [P23229](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15220

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD49f Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM