

Anti-CD53 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0273-FITC

Aliases

CD53; CD53 Molecule; TSPAN25; CD53 Antigen; MOX44; Leukocyte Surface Antigen CD53; Cell Surface Glycoprotein CD53; Tetraspanin-25; Tspan-25; Antigen MOX44 Identified By Monoclonal Antibody MRC-OX44; Transmembrane Glycoprotein; CD53 Tetraspan Antigen; Cell Surface Antigen; CD53 Glycoprotein

Background

Gene Name: CD53

NCBI Gene Entry: [963](#)

UniProt Entry: [P19397](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15230

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD53 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM