

Anti-CD54 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0275-FITC

Aliases

ICAM1; Intercellular Adhesion Molecule 1; CD54; BB2; Major Group Rhinovirus Receptor; ICAM-1; Intercellular Adhesion Molecule 1 (CD54), Human Rhinovirus Receptor; Epididymis Secretory Sperm Binding Protein; Cell Surface Glycoprotein P3.58; Human Rhinovirus Receptor; CD54 Antigen; P3.58

Background

Gene Name: ICAM1
NCBI Gene Entry: [3383](#)
UniProt Entry: [P05362](#)

Application Information

Clonality: Mouse monoclonal antibody
Clone ID: 24GB15235
Species Reactivity: Human
Applications Tested: Flow Cytometry (FCM)
Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD54 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM