

# Anti-CD57 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0281-FITC**

## Aliases

B3GAT1; Beta-1,3-Glucuronyltransferase 1; Galactosylgalactosylxylosylprotein 3-Beta-Glucuronosyltransferase 1; Glucuronosyltransferase P; GlcAT-P; UDP-GlcUA:Glycoprotein Beta-1,3-Glucuronyltransferase; EC 2.4.1.135; GLCATP; CD57; LEU7; CD57 Antigen; LEU7 Antigen; GlcUAT-P; GLCUATP; HNK1; NK-1; NK1

## Background

Gene Name: B3GAT1

NCBI Gene Entry: [27087](#)

UniProt Entry: [Q9P2W7](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14710

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgM

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD57 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)