

Anti-CD58 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0282-FITC

Aliases

CD58; CD58 Molecule; LFA3; CD58 Antigen, (Lymphocyte Function-Associated Antigen 3); Lymphocyte Function-Associated Antigen 3; Surface Glycoprotein LFA-3; Ag3; CD58 Antigen; LFA-3

Background

Gene Name: CD58

NCBI Gene Entry: [965](#)

UniProt Entry: [P19256](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14715

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD58 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM