

Anti-CD59 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0283-FITC

Aliases

CD59; CD59 Molecule (CD59 Blood Group); Membrane Attack Complex Inhibition Factor; Membrane Inhibitor Of Reactive Lysis; CD59 Glycoprotein; 16.3A5; P18-20; HRF20; MACIF; MIC11; MSK21; EJ16; EJ30; EL32; G344; MIN1; MIN2; MIN3; MIRL; CD59 Antigen P18-20 (Antigen Identified By Monoclonal Antibodies 16.3A5, EJ16, EJ30, EL32 And G344); CD59 Molecule, Complement Regulatory Protein; 20 KDa Homologous Restriction Factor; MEM43 Antigen; 1F5 Antigen; Protectin; HRF-20; MAC-IP; Surface Antigen Recognized By Monoclonal Antibody 16.3A5; Membrane Attack Complex (MAC) Inhibition Factor; CD59 Antigen, Complement Regulatory Protein; Lymphocytic Antigen CD59/MEM43; Human Leukocyte Antigen MIC11; T Cell-Activating Protein; CD59 Blood Group Antigen; MAC-Inhibitory Protein; Ly-6-Like Protein; CD59 Molecule; CD59 Antigen; MEM43; 1F5

Background

Gene Name: CD59

NCBI Gene Entry: [966](#)

UniProt Entry: [P13987](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15245

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD59 Mouse Monoclonal Antibody, FITC Conjugated


PAGE 2

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM