

Anti-CD62L Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0287-FITC

Aliases

SELL; Selectin L; LAM-1; PLNHR; CD62L; LYAM1; LSEL; LAM1; LNHR; Leukocyte-Endothelial Cell Adhesion Molecule 1; CD62 Antigen-Like Family Member L; Leukocyte Surface Antigen Leu-8; Lymphocyte Adhesion Molecule 1; Lymph Node Homing Receptor; L-Selectin; Gp90-MEL; Lyam-1; LECAM1; HLHRc; Leu-8; TQ1; Leukocyte Adhesion Molecule 1; Pln Homing Receptor; CD62L Antigen; LEU8

Background

Gene Name: SELL

NCBI Gene Entry: [6402](#)

UniProt Entry: [P14151](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15255

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD62L Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM