

# Anti-CD73 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0296-FITC**

## Aliases

NT5E; 5'-Nucleotidase Ecto; CALJA; CD73; ENT; NT5; EN; IMP-Specific 5'-Nucleotidase; Thymidylate 5'-Phosphatase; 5'-Deoxynucleotidase; Ecto-5'-Nucleotidase; 5'-Nucleotidase; EC 3.1.3.5; 5'-NT; NTE; 5'-Nucleotidase, Ecto (CD73); Purine 5-Prime-Nucleotidase; 5' Nucleotidase (CD73); CD73 Antigen; EC 3.1.3.35; EC 3.1.3.89; EC 3.1.3.91; EC 3.1.3.99; E5NT; NT

## Background

Gene Name: NT5E

NCBI Gene Entry: [4907](#)

UniProt Entry: [P21589](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14760

Species Reactivity: Human, african green, baboon

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD73 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)