

Anti-CD74 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0297-FITC

Aliases

CD74; CD74 Molecule; DHLAG; CLIP; CD74 Antigen (Invariant Polypeptide Of Major Histocompatibility Complex, Class II Antigen-Associated); CD74 Molecule, Major Histocompatibility Complex, Class II Invariant Chain; HLA Class II Histocompatibility Antigen Gamma Chain; Class II MHC-Associated Invariant Chain Peptide; HLA-DR Antigens-Associated Invariant Chain; Gamma Chain Of Class II Antigens; Ia-Associated Invariant Chain; MHC HLA-DR Gamma Chain; HLA-DR-Gamma; Ia Antigen-Associated Invariant Chain; CD74 Antigen; Ia-GAMMA; HLADG; P33; II; Ii

Background

Gene Name: CD74

NCBI Gene Entry: [972](#)

UniProt Entry: [P04233](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15270

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD74 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM