

# Anti-CD80 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0300-FITC**

## Aliases

CD80; CD80 Molecule; CD28LG1; CD28LG; B7.1; B7-1; CD80 Antigen (CD28 Antigen Ligand 1, B7-1 Antigen); T-Lymphocyte Activation Antigen CD80; B-Lymphocyte Activation Antigen B7; CTLA-4 Counter-Receptor B7.1; Activation B7-1 Antigen; LAB7; BB1; B7; Costimulatory Molecule Variant IgV-CD80; Costimulatory Factor CD80; CD80 Molecule; CD80 Antigen

## Background

Gene Name: CD80

NCBI Gene Entry: [941](#)

UniProt Entry: [P33681](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15275

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT


SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM

# Anti-CD80 Mouse Monoclonal Antibody, FITC Conjugated



## SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

## ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)